Materialography

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

Light curing device

Standard mode with 4 bulbs

# for polymerization of light-curing Kulzer resins

## ENGLISH

The Technotray POWER light curing device has been especially developed for polymerization of light-curing Kulzer resins (e.g. Technovit 2000 LC and Colorit).

The Technotray POWER can be equipped with up to 6 blue light tubes with 9 Watt each and is a high performance light-curing device that guarantees a rapid, intensive and homogeneous polymerisation.

The blue lights reliably cure layers of

- transparent Technovit resins with a thickness of up to 30 mm.
- Basic Colorit colours up to 0.8 mm and Deep Colorit colours with a layer thickness of up to 0.2 mm. With the exception of Colorit YELLOW, Technotray POWER reliably cures all Deep colours with a layer thickness of up to 0.8 mm and all Basic colours with a layer thickness of up to 0.2 mm.

The device's internal dimensions of WxDxH 170 x 160 x 120 mm allow for a spacious and homogeneously illuminated curing chamber, which can hold up to 4  $\emptyset = 50$  mm and up to 9  $\emptyset = 25$  mm embedding moulds [4].

### Starting your Technotray POWER for the first time

1. Check whether your mains voltage matches the required nominal voltage specified on the device's type plate. If necessary, set the device to the relevant mains voltage using the selector switch on the back [2].
2. Depending on the area of application, the device can be operated with either 4 or 6 light tubes (see "Curing" for more information).
  - Operation with 4 light tubes: 3 x lower sockets / 1 x upper sockets (middle) [3].
  - Operation with 6 light tubes: 3 x lower sockets / 3 x upper sockets [4].
  - a. Pull out the light curing tray as far as possible.[1]
  - b. Insert the light tubes into the light sockets until they engage and you can hear a click [3].
  - c. Close the light curing tray again.

3. Connect the device to the mains power. Power cord is included.
4. Press the power switch on the back of the device, which activates the lights and the timer in standby mode (the power switch will light up green when the device is switched on).
5. To open the curing chamber, pull out the light curing tray [1].
6. Select the curing time by pressing the timer button on the back of the device's housing (pressing the button 1 x = 5 min, pressing it 2 x = 10 min, pressing it 3 x = continuous operation) [5].
7. Fully push the curing tray back into the device again. On closing the device's door, it will automatically start the curing process for the selected time. The device automatically switches off if the tray is opened during operation. As soon as the tray is closed again, it will start the curing process anew for the selected time.

### Curing

1. Press the power switch on the back of the device, which activates the lights and the timer in standby mode (the power switch will light up green when the device is switched on).
2. To open the curing chamber, pull out the light curing tray [1].
3. Place the sample onto the glass base of the tray
4. Select the curing time by pressing the timer button on the back of the device's housing (pressing the button 1 x = 5 min, pressing it 2 x = 10 min, pressing it 3 x = continuous operation) [5].
5. Close the curing tray

# Technotray® POWER

## Light curing device

The Technotray POWER will automatically start the curing process for the selected time as soon as the curing tray has been closed. The lights will automatically switch off if the tray is opened while curing is in progress, and start up anew for the selected curing time as soon as the tray is closed again. The lights automatically switch off when the set curing time has elapsed.

### Curing examples Technovit 2000 LC

Standard Curing samples with covering varnish (Duration: 2 x 10 min) with the Technotray POWER

Place the sample inside a semi-transparent PE embedding mould. If required, the sample can be held in place with Technovit 2000 LC fixing paste.

Now fill the embedding mould with Technovit 2000 LC up to about 5 mm below the mould's top edge and place it on the Technotray POWER's curing tray. Set the timer to 10 min and close the tray. Once this curing cycle has elapsed, cover the sample with approx. 3 mm of Technovit 2000 LC covering varnish. This varnish will prevent a dispersion layer from forming on the bottom of the sample. Next, continue curing the sample for another 10 min.

Once the curing time has elapsed, remove the sample, leave to cool and then remove it from the mould.

Curing in layers (Duration: 5 min cycles) with the Technotray POWER

Embedding samples on a layer-by-layer basis will allow you to use significantly lower curing temperatures. Place the sample inside a semi-transparent PE embedding mould. If required, the sample can be held in place with Technovit 2000 LC fixing paste.

Now fill the embedding mould with approx. 2 - 3 mm Technovit 2000 LC, place it on the Technotray POWER's curing tray and cure for 5 min. Next, leave the sample to cool and then repeat this process until the sample is adequately embedded in Technovit 2000 LC.

The Technotray POWER can be used to simultaneously cure the following number of embedding moulds:

Ø 50 mm = 4 moulds; Ø 30 mm = 9 moulds  
Ø 40 mm = 5 moulds; Ø 25 mm = 9 moulds

### Curing temperatures

40 mm mould without sample completely filled:  
max. 120 °C

40 mm mould cured in individual layers (2 mm):  
5 min exposure: 70 – 90 °C

### Curing example for Colorit decorative resins

Standard curing process (Duration: 2 x 5 min) with the Technotray POWER device

Colorit decorative resins have to be cured using all 6 light tubes to ensure even polymerisation and homogeneous exposure to the curing light.

Place the coloured decorative item/items onto the curing tray. Set the timer to 5 min and close the tray. Once the curing time has elapsed, repeat the process. Once the decorative item/items have been cured for 2 x 5 min, they can be removed and processed as required.

# for polymerization of light-curing Kulzer resins

ENGLISH

## Changing the light tubes

We recommend changing the device's light tubes (OSRAM Dulux S9 W/71-blue) once a year or after 1000 operating hours or 20000 curing cycles.

Caution: Only change the light tubes when the device is cool!

To change the lights, perform the following actions in the order stated:

1. Switch off the device at the switch on the back.
2. Disconnect the mains plug.
3. Pull out the light curing tray as far as the lamps are completely accessible.
4. Next, start by removing the light tubes from the bottom row from their sockets, and then the light tubes from the top row.
5. Now insert the new light tubes into the sockets in reverse order.
6. Reconnect the mains power plug and switch the device on again on the rear of the unit.

## Cleaning the device

Always remove the mains plug before cleaning the device! The device must only be cleaned using a damp cloth. When doing so, always make sure that no fluids get into the device. Do not use any aggressive cleaning agents or solvents. Spare fuses can be obtained from specialist dealers.

## Technical data

- Mains voltage: 230 Volt, 50 Hz
- Output: max. approx. 60 watt
- 6 light tubes with 9 watt each
- Type of lamp: Blue light tubes

- Timer with 3 time settings: 5 minutes, 10 minutes and continuous operation
- Internal dimensions: WxDxH 170 x 160 x 120 mm
- Starts automatically when the tray is closed
- High-quality aluminium reflectors
- Light tube service life: approx. 1,000 operating hours or 20,000 cycles
- Housing colour: white/grey
- Dimensions: LxWxH 270 x 240 x 170 mm
- Weight: approx. 3 kg

Subject to technical changes.

## Package contents

Device and enclosed loosely [6]:

- 6 x light tubes
- Power cord
- Operating manual



## Disclaimer

We do not accept any liability for any use of the device other than specified in the operating manual or following unauthorised modifications.

Kulzer GmbH  
Division Technik  
Philipp-Reis-Straße 8/13  
61273 Wehrheim  
Tel. 06181-9689-2571 or 2574  
technik.wehrheim@kulzer-dental.com  
www.kulzer-technik.de

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